



ELISA ENZYME LINKED IMMUNOSORBENT ASSAY

Microwell Method

PROLACTIN

Cat. No. Z05303

P r o d u c t I n s e r t



Enzyme Linked Immunosorbent Assay for the **quantitative**
determination of Prolactin in human serum or plasma.

Microwell Method - 96 wells

(12 x 8-well Antibody coated Strips)

Individual breakaway

GENERAL INFORMATION

☐ **Wavelength**

Measurement Filter: 450/620 nm

☐ **Incubation Time**

80 minutes at RT (60/20)

☐ **Enzyme Conjugate**

HRP (Horseradish Peroxidase)

☐ **Substrate**

TMB (3,3',5,5'-Tetramethyl-benzidine)

☐ **Sample**

Serum or Plasma

☐ **Stability of Samples**

undiluted: 2 days at 2-8°C; for longer storage at - 20 °C

☐ **Calibration Range**

0 – 200 ng/ml

☐ **Sensitivity**

1 ng/ml

☐ **Shelf life and Stability of Kit Components**

Kit: 18 months from production date.

Kit Components: see expiration date on the label

KIT COMPONENTS

Microwell plate	12x	8 well strips with breakaway microwells coated with anti-monoclonal PRL.
Calibrators	5	vials, A 1 ml, B-E 0,15 ml, ready to use . Aprox. 0-5-20-75-200 ng/ml , see labels
Enzyme Conjugate	1	vial of anti-monoclonal PRL-HRP conjugate, 11 ml, ready to use .
Substrate Solution	1	vial of H ₂ O ₂ , TMB 0.25 g/l, 11 ml, ready to use .
Stop Solution	1	vial of sulphuric acid 0.15 mol/l, 11 ml, ready to use .

MATERIALS REQUIRED BUT NOT PROVIDED

- ☐ Deionized or distilled water
- ☐ Graduated cylinders and beakers
- ☐ Wash trays
- ☐ **Macropipettes** capable of delivering 5 µl to 1000 µl. We recommend Socorex Digital Macropipette.
- ☐ **Multichannel Micropipette**; we recommend Socorex Digital Multichannel Micropipette.
- ☐ **Stepper**; we recommend Socorex Stepper.
- ☐ **Microplate reader** capable of reading absorbance values at 450 nm. If dual wavelength microplate reader is available, the reference filter should be set at 600-690 nm.
- ☐ **Automatic microplate washer** capable of dispensing 200 µl.

Microplate reader and microplate washers are available from Dialab Company.

SUMMARY AND EXPLANATION

Prolactin (PRL) is a polypeptide hormone secreted by the anterior pituitary gland. Release of prolactin is pulsatory, and regulated by hypothalamic releasing and inhibitory factors. High prolactin levels are observed during pregnancy and the postpartum-lactation period. Other physiological states associated with high prolactin levels include physical and emotional stress, sleep and hypoglycaemia. Hyperprolactinaemia can cause suppression of gonadal function. Prolactin measurement therefore forms essential part of investigations of infertility. In woman, symptoms that may accompany hyperprolactinaemia include nonpuerperal galactorrhoea, amenorrhoea and other menstrual disorders. In men it may be associated with loss of libido and impotence. Pathological causes of hyperprolactinaemia include prolactin-secreting macro- and microadenomas, hypothyroidism and renal failure.

TEST PRINCIPLES

The ELISA test is performed as an indirect solid phase **sandwich-type** immunoassay. Microwells are coated with anti-monoclonal PRL followed by blocking the unreacted sites to reduce non-specific binding.

Step 1 Prolactin Antigens present in calibrators and patient samples bind to the coated antibody.

Step 2 The Antigen-Antibody complex is reacted with enzyme (HRP) labeled anti-monoclonal PRL conjugate resulting in the PRL antigen being sandwiched between the solid phase antibody and the enzyme conjugate.

Step 3 The enzyme converts added substrate (TMB) to form a colored solution.

Step 4 The intensity of color change, which is proportional to the concentration of Antibodies present in the samples is read by a microplate-reader at 450 nm. Results are expressed in ng/ml.

EXPECTED VALUES

Each laboratory must establish its own normal ranges based on patient population. The serum PRL values are comprised in the following intervals:

Sample	Range ng/mL
Children 0-14 years	2,1 – 17.0
Male:	2.1 – 17.7
Female:	
14 – 55 years	2.8 - 29.2
pregnant	9.7 – 208.5
menopause	1,8 – 20.3

Some of the female population tested in this group were probably using oral contraceptives, which may affect results.

REAGENTS

Storage

☐ Store all reagents at 2° - 8°C. Do not freeze!

Precautions

- ☐ Instructions should be followed exactly as they appear in this kit insert to ensure valid results.
- ☐ Avoid contact with the **TMB (3,3',5,5'-Tetramethyl-benzidine)**. If TMB comes into contact with skin wash thoroughly with water and soap.
- ☐ The stop solution contains **sulphuric acid**. If it comes into contact with skin, wash thoroughly with water and seek medical attention.
- ☐ Avoid contact between the buffered **peroxide** solution and easily oxidized materials; extreme temperatures may initiate spontaneous combustion.
- ☐ Do not use beyond expiration date on the label.
- ☐ Do not use if reagent is not clear or if a precipitate is present.
- ☐ Do not interchange kit components with those from other sources other than the same catalog number from DIALAB.
- ☐ Follow good laboratory practices to minimize microbial and cross contamination of reagents when handling.
- ☐ All human derived components used have been tested for HBsAg, HCV, HIV-1 and 2 and HTLV-I and found negative by FDA required tests. However, human blood derivatives and patient specimens should be considered potentially infectious. Follow good laboratory practices in storing, dispensing and disposing of these materials.

SPECIMEN COLLECTION AND HANDLING

- ☐ Only **Serum or Plasma** specimens should be used in this procedure. The patients need not to be fasting, and no special preparations are necessary.
- ☐ Grossly hemolyzed, lipemic or microbially contaminated specimens may interfere with the performance of the test and should not be used. Neither Bilirubin or Hemolysis have significant effect on the procedure.
- ☐ Store specimens at **2°- 8°C** for up to a maximum of **2 days**. For longer storage, specimens should be frozen. Avoid repeated freezing and thawing of samples.

PROCEDURE

Procedural Notes

- ☐ Before starting with the assay read carefully the product insert.
- ☐ Let specimens and test reagents equilibrate at room temperature before starting with the test procedure. Return all unused specimens and reagents to refrigerator immediately after use.
- ☐ Remove required microwell strips from the pouch and carefully reseal the pouch to prevent condensation in the unused wells. Return pouch immediately to refrigerator.
- ☐ Good washing technique is critical. For manual washing, fill each microwell with 200 µl distilled water. Discard the fluid by inverting and tapping out the contents of each well or by aspirating the liquid from each well. To blot at the end of the last wash, invert strips and tap the wells vigorously on absorbent paper towels. For automatic washers, program the washer as per manufacturer's instructions.
- ☐ Use a multichannel pipette capable of delivering 8 wells simultaneously. This speeds the process and provides for a more uniform incubation time.
- ☐ For all steps, careful control of timing is important. The start of all incubation periods begins with the completion of reagent addition.

Preparation of sample

Usually no dilution necessary; dilute samples with concentrations above 200 ng/ml 1:1 with Calibrator A.

Test Procedure

Step 1 Label protocol sheet to indicate sample placement in the wells according to the following figure. **5 calibrators** (SA-SE) and **1 Blank** should be included. The user has the option to run Patient Samples (P) in duplicate.

	1	2	3	4	5	6	7	8	9	10	11	12
a	B	SD	P3									
b	SA	SE	P4									
c	SA	SE	P4									
d	SB	P1	P..									
e	SB	P1	P..									
f	SC	P2										
g	SC	P2										
h	SD	P3										

Calibrator	Conc. ng/ml
SA	0
SB	5
SC	20
SD	75
SE	200

Step 4 Remove the required microwells from pouch and return unused strips in the sealed pouch to refrigerator. Securely place the microwells into the extra provided holder .

Step 5 Pipette **10 µl** of **Calibrators and Patient Samples** into the wells.

Step 6 Add **100 µl** of diluted **Enzyme Conjugate** to the wells except for Blank well and incubate **60 minutes** at room temperature.

Step 7 Discard the contents of the microwells and wash the wells with **200 µl** distilled water. Repeat 2 x the washing procedure by draining the water completely.

Step 8 Pipette **100 µl of Substrate Solution** into each microwell in the same order and timing as for the Enzyme Conjugate, Blank well included.

Step 9 Incubate **20 minutes** at room temperature in the dark.

Step 10 Add **100 µl of Stop Solution** into each microwell using the same order and timing as for the addition of the Substrate Solution.

Step 11 Read absorbance of each microwell at **450 nm(reference 620 nm)** **against blank** using a microplate reader. The developed color is stable for at least 30 minutes. Read optical densities during this time.

TEST EVALUATION

Mean absorbance and relative percentage

1. Calculate the mean of the absorbances (E_m) corresponding to the single points to the standard curve and of each sample.
2. Subtract the mean absorbance value of the zero calibrator from the mean absorbance values of calibrators and samples.
3. Draw the standard curve on log-log graph paper by plotting absorbance values of standard against appropriate PRL concentration.
4. Read off the PRL concentrations of the control and samples.

LIMITATIONS OF THE PROCEDURE

The assay should not be performed on grossly hemolyzed, microbially contaminated or lipemic samples. This method should be used for testing human serum samples only.

Q.C. PARAMETERS

Maximum Absorbance (calibrator E) = OD (200 ng/ml) = min. 50% of the value stated in quality certificate.

PERFORMANCE CHARACTERISTICS

Sensitivity

The minimal detectable concentration of Human Prolactin by this assay is estimated to be 1.0 ng/ml.

Specificity

The cross reaction of the antibody calculated at 50% according to Abraham are shown in the table:

hPRL	100.0 %
hLH	< 1 %
hCG	< 1 %
hFSH	< 1 %

Precision

a. Intra Assay variation

Within-run precision was determined by replicate determination of three different control in one assay. The within assay variability is shown below:

Sample	1	2	3
Number of replicates	12	12	12
Mean PRL (ng/ml)	6,7	16,4	36,5
Coef. of Variation (%)	6,2	8,5	5,25

b. Inter Assay variation

Between-run precision was determined by replicate determination of three different controls in one assay. The between assay variability is shown below:

Sample	1	2	3
Number of replicates	12	12	12
Mean PRL (mIU/ml)	6,7	15,52	36,43
Coef. of Variation (%)	6,4	8,6	3,9

7.4 Recovery

Expected conc.	Observed conc.	Recovery
9,67	9,83	98,4
34,94	35,75	97,8

7.5 Linearity

Two patient samples were serially diluted with zero standard in a linearity study. The average recovery was 102.2 %.

Dilution	Exp. Conc	Obs. Conc	Recovery
		32,5	
Dil. 1 / 2	16,65	16,25	102,5
Dil. 1 / 4	7,84	8,12	96,6
Dil. 1 / 8	4,36	4,06	107,4

REFERENCES

1. Wisdom, G.B. Clin. Chem. 22/8 1243 - 1255 (1976)
2. Shome, B. and Parlow, A.F., J. Clin. Endocr. Metab.,39, 199 – 202 (1974).
3. Shome, B. and Parlow, A.F., J. Clin. Endocr. Metab.,39, 203 – 205 (1974).
4. Uotila, M. Ruoslahti, E. Engvall, E., J. Immunol. Methds, 42, 11-15 (1981)
5. Acosta, A.A.M.D. and Wright, G.L., Journal of Clinical Immunoassays, 6, 41 (1983).
6. Cohen, K. L., Metabolism, 26 1165-1177 (1977)

ELISA Enzyme Linked Immunosorbent Assay



DIALAB Ges. m.b.H.

Produktion und Vertrieb von chemisch – technischen Produkten und Laborinstrumenten Gesellschaft m.b.H.
A-2351 Wiener Neudorf – IZ NÖ Süd, Hondastrasse, Objekt M55 - Phone + 43 (0) 2236660910-0
Fax + 43 (0) 2236660910-30

E – Mail: office@dialab.at URL: www.dialab.at